

TRACHEOSTOMY REVIEW AND MANAGEMENT SERVICE CLINICAL PROCEDURE

PASSY MUIR VALVE (PMV) USE IN SPONTANEOUSLY BREATHING PATIENTS

Staff this document applies to:

Medical Staff, Nurses, Speech Pathologists, Physiotherapists on all campuses, including ICU and in the community.

State any related Austin Health policies, procedures or guidelines:

[Passy-Muir Valve \(PMV\) Use In Line with the Ventilator](#)

[Tracheostomy Cuff Deflation and Reinflation](#)

[Tracheostomy Clinical Procedure - Decannulation](#)

[Clinical Instruction Sheet - Use of Portex Suctionaid Tracheostomy Tube](#)

[Tracheostomy Learning Package - Use of the Passy-Muir Valve](#)

[Stridor in Adults](#)

Definition:

The Passy Muir Valve (PMV) is a one-way valve that opens upon inspiration and closes completely upon expiration. It fits the 15mm hub of any standard tracheostomy tube. A valve is always used with the cuff fully deflated.

Rationale:

- A one-way valve restores airflow through the upper airway which facilitates voicing, coughing, swallowing, return of sensation, and smell.
- Using a one-way is an important step in the decannulation process.

Clinical Alert:

The cuff must always be deflated when a Passy Muir Valve (PMV) is used. Failure to deflate the cuff can result in immediate respiratory compromise/distress and/or death.

- The PMV should be removed when:
 - The patient is sleeping
 - The patient exhibits respiratory distress +/- stridor including increased work of breathing
 - Patient appears distressed or is in discomfort specifically when the PMV is on
 - The patient becomes unable to produce voice or there is a sudden change in their voice
- Do not place the PMV under the following circumstances:
 - Increased or copious secretions
 - Sub-optimal humidification
 - Upper airway obstruction
 - Difficulty passing a suction catheter
 - Severe coughing

- Unstable cardiorespiratory status: The specific cardiorespiratory parameters that are acceptable may vary depending on the patient population. Liaise with the treating medical team regarding the acceptable cardiorespiratory parameters.
- With a Bivona Fome-Cuf tracheostomy tube
- If a patient arrives with a speaking valve from another manufacturer, remove the valve and contact the treating Speech Pathologist to assess the patient with a PMV
- Only staff trained in placement of PMV may use PMV. Completing the PMV e-learning package on ATLAS is highly recommended.

Who is authorised to perform this procedure?:

- Following consultation with the treating Medical staff, the Speech Pathologist will conduct assessment for initial PMV placement with assistance from Physiotherapist and/or Specialist Nurse
- If outside business hours, senior medical staff in ICU with training in valve placement can place the PMV for patient assessment. During business hours the Speech Pathologist will conduct formal assessment and establish the schedule for routine use if indicated

Expected Outcome:

- The patient will safely and comfortably wear the PMV
- The patient will be able to voice, cough, clear secretions and swallow saliva with increased ease
- The patient will routinely use the PMV according to the schedule established by the Speech Pathologist

Equipment:

- Routine tracheostomy personal protective equipment (PPE)
 - Clean gloves
 - Safety shield, goggles or glasses
 - Disposable apron
 - Surgical mask
- Functioning suctioning equipment, suction catheters and Yankuer
- 10 ml syringe
- PMV (# 007 Aqua is the valve model most frequently used at Austin Health)
- Pulse oximetry during initial sessions or if patient status changes

Procedure:

- Explain procedure to patient
- Note baseline measures: HR, breathing patterns, RR and SpO2 if using oximeter
- If a Suctionaid tracheostomy tube is in situ, use a syringe to aspirate the Suctionaid line and mouth to clear above-cuff secretions
- Suction via tracheostomy
- Completely deflate cuff and suction via tracheostomy
- Check pilot balloon to ensure cuff is completely deflated
- Secure the tracheostomy and place PMV on the tracheostomy with quarter clockwise turn
- If patient's voice sounds wet or gurgly there may be secretions in the upper airway. Encourage the patient to bring those to the mouth

- Clear mouth with Yankauer (if required)
- If patient cannot voice or demonstrates changes to heart rate, breathing patterns, respiratory patterns or SpO₂, or is visibly distressed, do not use the PMV. Contact the Speech Pathologist

Post Procedure:

- Ensure patient is comfortable and able to voice as expected
- Monitor patient closely for intolerance of the valve including increased work of breathing, anxiety, increased coughing
- Remove the valve as per documented schedule or at any time if the patient does not tolerate the valve. Document duration and tolerance with PMV use in the medical record

Author/Contributors:

Tanis Cameron (TRAMS Manager, Senior Speech Pathologist)

Prue Gregson, Adrianna Griffith & Charissa Zaga (Senior Speech Pathologists)

Tracheostomy Policy and Procedure Committee

Updated October 2020

Legislation/References/Supporting Documents:

Kazandjian M.S., Dikeman, KJ. (2008). Communication options for tracheostomy and ventilator-dependent patients in *Tracheotomy: Airway Management, Communication and Swallowing*, 2nd ed. Edited by Eugene N. Myers and Jonas T. Johnson Chapter 11 p 187-214 Plural Publishing, Inc., San Diego, CA.

Morris, L. (2010). *Tracheostomies. The Complete Guide*. Springer Publishing Company. New York.

Morris, L. L., Bedon, A. M., McIntosh, E., & Whitmer, A. (2015). Restoring speech to tracheostomy patients. *Critical care nurse*, 35(6), 13-28.

Passy- Muir. (2012). *Passy- Muir tracheostomy & ventilator swallowing and speaking valves instruction booklet*. CA: Passy-Muir Inc.

Ross, J., White, MF. (2003) Removal of the tracheostomy tube in the aspirating spinal cord injured patient. *Spinal Cord* 41:636-642

Russell, C., Matta, B. (2004). *Tracheostomy: A multiprofessional handbook*. London: Greenwich Medical Media Limited.

Speed, L., & Harding, KE. (2013). Tracheostomy teams reduce total tracheostomy time and increase speaking valve use: A systematic review and meta- analysis. *Critical Care*. 28(2): 216.e1-10.

Suiter, DM., McCullough, H., Powell, W. (2003). Effects of cuff deflation and one-way tracheostomy speaking valve placement on swallow physiology *Dysphagia* 18: 284 -292.

Suiter, DM., Leder, SB. (2007). Contribution of tracheotomy tubes and one-way valves to swallowing success. *Topics in Geriatric Rehabilitation* 23(4): 341-351.

Dikeman K, Kazandjian M. *Communication and swallowing management of tracheostomized and ventilator-dependent adults*. 2nd ed: Delmar Learning; 2003.

Hess, D. (2005). Facilitating speech in the patient with a tracheostomy. *Resp Care* 50(4) 519-525

Tippett, D. (2000) *Tracheostomy and ventilator dependency; management of breathing, speaking and swallowing*: Thieme Medical Publishers.

Authorised/Endorsed by:

A/Prof Stephen Warrillow, A/Prof Mark Howard, Joanne Sweeney, A/Prof Sue Berney
Tracheostomy Policy and Procedure Committee
Clinical Nursing Standards Committee

Primary Person/Department Responsible for Document:

TRAMS